



## ELECTRO SURGICAL UNIT



New Micro controller based electro surgical units are available with outstanding performance of monopolar (cut & coag) & Bipolar.

In many types of surgery where high precision & reliability are essential, microcontroller control provides the safety & performance during surgery.

These units have pure cut are exceptionally smooth & start Promptly even in irrigated procedure.

The blend gives a good hemostasis effect

than cutting & spray coagulation which rivals the best spark gap generators & with minimum cutting effect. The isolated bipolar output, which has a non sparking characteristic, is ideal for microsurgery, neurosurgery, laparoscopy & other application.

DISPLAY	: 7 SEGMENT LED
DIMENSION (WXHXD)	: 290 MM X 105 MM X 260 MM
WEIGHT	: 4 KG.
KEYBOARD	: FEATHERTOUGH
INPUT SUPPLY VOLTAGE	: 230VAC + 15% @ 50/60HZ

