




# TROCAR SYSTEMS



Cat #	Description
TRO-20-F	Trocar Sleeve with Trocar Pyramid Tip, 20mm, Flap
TRO-10-F	Trocar Sleeve with Trocar Pyramid Tip, 11mm, Flap
TRO-10-M	Trocar Sleeve with Trocar Pyramid Tip, 10mm, Multifunction
SS-10P	Spiral Sleeve compatible with TRO-10-MP and TRO-10-FP
HA-10	Hassan Cone compatible with TRO-10M and TRO-10-F
SET-10P	Safety Trocar compatible with TRO-10-MP and TRO-10-FP
BT-10P	Trocar Blunt Tip compatible with TRO-10-MP and TRO-10-FP

Cat #	Description
TRO-07-F	Trocar Sleeve with Trocar Pyramid Tip, 7mm, Flap
TRO-05-F	Trocar Sleeve with Trocar Pyramid Tip, 5.5mm Flap
TRO-05-FP	Trocar Sleeve Light Body with Trocar Pyramid Tip, 5.5mm, Flap
TRO-05-M	Trocar Sleeve with Trocar Pyramid Tip, 5.5mm, Multifunction
SS-05P	Spiral Sleeve compatible with TRO-05-FP
TRO-04-M	Trocar Sleeve with Trocar Pyramid Tip, 3.5mm, Multifunction

Note: All above Trocar Sleeves / Trocar are available in extra long length on request.

## REDUCERS



Cat #	Description
CRS-20-10	Reducer Sleeve 20-10
CRS-10-05	Reducer Sleeve 11-5
CRS-10-15	Clip on Reducer 11-5 fitting TRO-10-F & TRO-10-M
CRS-07-05	Reducer Sleeve 7-5
CRS-05-03	Reducer Sleeve 5-3

## VERESS NEEDLE



Cat #	Description	Length
CVN-90	Veress Needle	90mm
CVN-100	Veress Needle	100mm
CVN-120	Veress Needle	120mm
CVN-150	Veress Needle	150mm
CVN-170	Veress Needle	170mm

## Instruments

- Handle made from plastic-Lightweight & Autoclavable
- Easy 3 Piece dismantling for better cleaning and sterilization
- Rotating & Dismantable.
- All handles, Insulated Shafts and Jaw Inserts are interchangeable.